# LIGAMENTOTAXOR

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### **DESCRIPTION**

This is a dynamic external digital fixator. It is made up of two 1.2 mm stainless steel pins, 4 plastic rods and two 0.8 mm diameter springs. To guide the double spring system, the surgeon insert a third 1.5 mm diameter pin. This third U-shaped pin interlocks the two distraction mechanisms.

For optimal protection, LIGAMENTOTAXOR also features 4 rod-protection tips and a single use drill guide to ensure the parallelism of the device's two 1.2 pins. For ideal application, LIGAMENTOTAXOR is combined with the 1.6TB MANOTTE pliers with which the pins can be bent at 90° (without mechanical stress to the bone). The pliers are also used to cut off the pins cleanly and to crimp the pin protection tips.

The device remains in place for up to a maximum of 60 days.

### **STERILIZATION**

The LIGAMENTOTAXOR contained in this package is sterile. It is sterilized by ionising radiation at a minimum dose of 25 kGy.

It can not be resterilized.

### **INDICATION**

LIGAMENTOTAXOR is a dynamic external fixation system intended for use in the treatment of complex fracture dislocations, stable or unstable dislocations, fracture luxations, and pilon fractures in the PIP joint. It is intended for use in treating PIP joint flexion contractures with an associated loss of range of motion. It helps restore proper digit alignment and range of motion while permitting early mobilization and normal anatomical movement.

For children, surgeon will determine if LIGAMENTOTAXOR can be used taking account his age and morphology as well as his ability to follow post-operative cares.

### **CONTRAINDICATION**

- Acute or chronic, local or systemic infections ; Fever.
- Known allergy to chromium and other components contained in stainless steel.
- Insufficient tissular cover.
- Infectious disease ; systemic and metabolic disorders.
- Drug or alcohol dependency.
- Severe osteoporosis, local bones tumors, severe stiffness of the joints.
- Poor blood supply of the soft tissues, insufficient bone stock.
- Evident patient intellectual incapacity to follow the surgeon's recommendations: non cooperative patient or with mental illness.
- Intense physical activity.

### **WARNING**

The LIGAMENTOTAXOR cannot be used with MRI as it may cause skin burns.

Patients must be fully informed concerning the objectives and modus operandi of distraction so that the procedure produces optimal results. Patients must be informed of warning signs which should lead them to call their surgeon for remedial action to be taken: appliance failure (loosening), traction failure, breakage, unexpected pain, infection, pus discharge.

### **SURGICAL PROCEDURE**

Implementation of appropriate surgical procedures and techniques are the sole responsibility of the medical profession.

## **PATIENT INFORMATION**

Some degree of risk is inherent to all surgical procedures and complications are sometimes unavoidable. It is the surgeon's responsibility to explain fully the benefit-to-risk ratio of the procedure under consideration. Potential complications and caveats must be discussed with patients before a decision is taken to proceed with surgery. The possibility of discontinuing treatment at any point must also be made clear to patients.

## **POSSIBLE COMPLICATIONS AND WARNINGS**

Complications which may arise as a result of all surgical procedures must be discussed with patients (infection, haematoma, accumulation of fluids, loss of sensitivity, pain, stiffness, finger deviation-deformity, Reflex Sympathetic Dystrophy Syndrome, subluxation, possible allergic reaction to medications used during surgery, nerve damage, patient intolerance to a foreign body, imperfect healing).

## **SINGLE USE**

Devices are designed for single use. Those removed after use must not be used again since re-cleaning and re-sterilization procedures may fail to eliminate with certainty biological waste such as blood, tissues and other materials that could contain resistant pathogenic germs.

## **PACKAGING**

The devices are supplied sterile, processed by strictly validated and controlled sterilization cycles. Sterility is verified to comply with current standards. Sterility of the device is only guaranteed if packaging is inviolate and intact. Removable labels to ensure traceability of the product are supplied and attached to primary packaging.

## **DIRECTIONS FOR USE**

Correct implementation of surgical procedures and techniques are at the discretion of members of the medical profession. Surgeons must individually evaluate the suitability of the procedure on the basis of currently accepted techniques, relying on their individual appreciation and experience. The following is a description of placement of the device as applicable to a PIP joint. Once the principle is understood, its applicability can easily be extended to other sites. An image intensifier must be made available.

The first KW is introduced transversally into the epiphyseal axis of P1, ideally central to the rotational axis of the joint. The second KW is introduced using the drill guide at P2 shaft level. Using the MANOTTE 1.6TB, the distal KW is bent at a 90° angle on both sides of the finger. This KW is then shortened as shown on the illustrated instructions for use. Slide the two plastic rods over the proximal KW, one on each side of the finger. Depending on the length of the patient's fingers, use either the long green plastic rods or the short blue ones. The long green rods are for long fingers. Furthermore they will be helpful for lateral X-rays.

The spring is then 'screwed' onto the section of the distal KW which is bent at 90° as shown on the photos. The spring is then pushed halfway forward. (The joint's reaction to ligamentotaxis can be tested by pulling the springs). The plastic rod is then curved slightly as it is introduced into the spring. The spring screwing procedure can then be completed in order to achieve the required distraction. N.B. Recommendation to facilitate the screwing in of the spring: either pull on the spring while it is screwed in, or alternatively, pull on one of the springs while the other spring is screwed in, and continue likewise. The proximal KW is bent at a 90° angle and cut on both sides of the finger, leaving sufficient space for the rod to rotate without rubbing against the skin. Finally, KW protective caps may be placed at each end, but this is not mandatory since the MANOTTE 1.6TB pliers cut cleanly and precisely. To stabilise the springs, a frame is made using 1.5 mm KW. The frame must be wide enough to allow each KW extrem-

ity to enter into one of the two springs whilst keeping them on the outside of the finger. To maintain the required degree of distraction, use two strips of adhesive tape covering the springs and the U-shaped KW. Finally make a note on the tapes of the number of turns of the springs between the two KWs once optimal distraction is achieved.

The degree of distraction depends on the clinical status and the condition of the joint. It may be necessary to apply considerable force on one side and less on the other. Only the surgeon, aided by the image intensifier, can decide on what degree of traction is required. To avoid dorsal subluxation at the base of P2, a third more proximal KW may be applied.

### **TO OPEN A PACKAGED STERILE PRODUCT**

- 1 Remove the outer packaging in an aseptic and sterile environment, letting the sealed internal packaging fall away gently.
- 2 Fix the registering portion of the internal label to the patient's case file.
- 3 Open the inner packaging in asepsis.

### **EXAMINATION AND HANDLING OF THE PRODUCT**

The product should be examined visually to detect any sign of damage.

### **CONTAMINATION**

AREX<sup>®</sup> cannot be held responsible for products contaminated by substances once the product has left its premises. Surface contamination (talcum powders, dust, dressings, wadding, oil) on the devices may trigger reactions since they are foreign bodies. Handle with care using surgical gloves (rinsed out to eliminate talcum powder) taking care to preserve asepsis at all times.

### **LIMITED WARRANTY, LIMITATION OF LIABILITY AND WAIVER OF OTHER GUARANTEES.**

The manufacturer guarantees that all necessary precautions were taken during manufacture of the products and undertakes to replace any product which, following the manufacturer's investigation, proves to have been defective when it was shipped. Patient selection, surgical procedures, treatment, post-operative stress and handling of the appliances are entirely the surgeon's responsibility. AREX<sup>®</sup> has no control over conditions of use nor on the good or bad effects of the product's use. AREX<sup>®</sup> will not be held responsible for loss or damages as a direct or indirect result of using the appliance. Destruction of waste will comply with current legislation.

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# INSTRUCTIONS FOR USE of LIGAMENTOTAXOR

For optimal use MANOTTE 1.6TB  
and Drill guide are mandatory.

QR-codes for video instructions :

LIGAMENTOTAXOR


<http://www.arex.fr/QRc.php?c=2>


MANOTTE 1.6TB

<http://www.arex.fr/QRc.php?c=8>




First, put in place proximal 1.2 mm KW [1]. Then use disposable drill guide [2] to ensure parallel placement of 2nd (distal) KW [3].





Open MANOTTE 1.6TB pliers as wide as possible and place under the hand. Bend the distal KW at a 90° angle on both sides of the finger [4].


Cut bent distal KW and place outwards along the finger. See photo above.

Slide blue plastic rods over proximal KW on both sides of the finger. Select the end of the spring with the **longest straight end**. Screw each spring completely over both distal KWs and blue plastic rods as shown *[5]*. The spring end will lock automatically. If extra traction is required, unscrew the spring away from the blue plastic rod. Pull distally while re-screwing.

Bend proximal 1.2 mm KW at 90° on each side [6]. Be sure to leave sufficient space to allow blue plastic rods to rotate without rubbing against the finger. Distribute the space left by MANOTTE on each side equally, as shown on photos 7, 8, 9.

Cut proximal KW as shown, then place KW protective caps as shown and finally crimp them into place.

To stabilise both springs, construct a “frame” using the 1.5 mm KW supplied [0.059in.] creating a U-shape large enough to maintain the springs away from the finger. If necessary, springs can be shortened before constructing the U-shaped frame.

Two strips of tape can be placed around both the springs and the U-shaped frame to secure selected distraction. The number of spring turns between the proximal and distal KWs can be noted on the tapes.

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