TCU

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DESCRIPTION

UNIVERSAL TENDON SPACER with an oval section for temporary implantation for use in the two stage reconstruction of the flexor and extensor tendons of the hand.

The radio-opacity of the product allows for X-ray viewing. Second stage surgery (material ablation and reconstruction of the tendons) is normally carried out between two and six months after implantation. The second stage operation is carried out after an appropriate pseudo-sheath for the tendon has developed around the implanted device.

Every package contains a single implant in sterile packaging. The product is for single use only. With its unique concept - progressively increasing dimensions, a diameter which increases from 3 mm to 6 mm and from 1.5 mm to 3 mm (smaller dimension / oval section) which reproduces the characteristic dimensions of the digital canal - this product is available in a single size which adapts naturally to the varying dimensions of the digital canal.

STERILITY

The UTS contained in this packaging is sterile. Sterilization is carried out using ionising radiation at a dose of at least 25 kGy. The product may not be re-sterilized.

INDICATIONS

Adherent or mutilated tendons following trauma or a primary repair which has failed. Absence of the tendon sheath. Non-functioning, mutilated or adherent tendon pulleys. Ruptured tendons.

The patient must be in good health with sufficient tissue cover. Sufficiently healthy joints, bone marrow and adequate vascular and neuro-muscular conditions are also necessary to facilitate a definite functional improvement using this reconstruction technique. The patient must be cooperative, especially concerning the need for rehabilitation measures.

CONTRA-INDICATIONS

Insufficient tissue cover. Infection. Severe stiffness of the joints. Poor nutrition of the soft tissues.

Associated injury to the extensor mechanism, nerve injury, involvement of several fingers, insufficient blood flow to the tissues.

WARNING

Cases with multiple prior surgical operations which have failed and with stiffness of the fingers are unlikely to obtain much benefit from this procedure. The patient must be regularly monitored between the two operating stages, both clinically and using X-rays, to follow the patient's progress and to anticipate complications. To avoid complications which may lead to a disinsertion of the distal attachment, synovitis in the pseudo-sheath of the tendon or exposure of the UTS by erosion, the second operation should not be delayed without good reason.

SURGICAL PROCEDURE

Full responsibility for carrying out the appropriate surgical procedures lies with the surgeon and medical staff.

PATIENT INFORMATION

All surgical procedures carry risks and can be accompanied by complications. The surgeon must clearly explain the benefits and risks accompanying the planned intervention. Possible complications and dangers must be discussed with the patient prior to making the decision to proceed with surgery. The possibility of having to remove the UTS at any moment must also be raised with the patient. The surgeon must inform the patient before the operation of potential complications and risks to which the patient will be exposed.

POSSIBLE COMPLICATIONS AND RISKS

Complications which may occur in the course of all surgical procedures must be discussed with the patient (infection, haematoma, accumulation of liquids, loss of sensation, possible allergic reaction to drugs used during the operation, lesion of nerves, intolerance by the patient of a foreign body, scarring).

Possible known complications that may arise when utilising the device:

- Secondary capsular contraction

The formation of a fibrous tissue capsule around the implant is a normal physiological response to a foreign body. The contraction of fibrous capsular tissue around the implant can lead to hardening, a sensation of discomfort or even pain. It has been suggested that haematomas which may occur after a surgical operation can increase the probability of a secondary capsular contraction.

- Rupture / breakage

Tear of distal attachment sites, pulley rupture, tenolysis/ additional surgery, implant breakage.

- Exteriorization

Insufficient tissue coverage can result in the exposure and exteriorization of the implant. A cutaneous rupture has been reported following the use of steroid drugs.

- Deterioration of the implant

Arrangements should be made for a supplementary implant to be available in the event of deterioration of the implant and/or accidental contamination during its implantation. Extreme care must be taken during the use and handling of the implant to minimise these eventualities.

CONTAMINATION

Care must be taken that surface contaminants such as talc, dust and oils from the skin do not come into contact with the implant.

SINGLE USE

These products are for single use only. Products which have been removed must not be reused, as the cleaning and re-sterilization procedures cannot be relied upon to remove all biological residues such as blood, tissue and other materials which could contain resistant pathogens.

INCORRECT SIZE, ASYMMETRY, BADLY PLACED SCAR TISSUE, BADLY POSITIONED IMPLANTS

These complications are often associated with surgical technique. The use of UTS should minimize such complications which are usually related to the choice of a tendon spacer which does not have the appropriate dimensions.

PACKAGING

The implants supplied are sterile and have been treated using strictly validated and controlled sterilization cycles. Sterility is controlled in accordance with current standards. The sterility of the implant is only guaranteed if its packaging remains intact and has not been damaged. Detachable labels to ensure traceability are supplied and are fixed to the outer packaging.

METHOD OF USE

The choice of surgical procedure and technique lies with the surgeon. Each surgeon must evaluate the suitability of the procedure based on currently accepted medical techniques, on an evaluation of the individual case and the surgeon's own experience. The use of UTS for a two-stage reconstruction of flexor and/or extensor tendons is the responsibility of the surgeon.

The surgeon must nevertheless ensure that the UTS occupies the maximum space allowed by the pulleys while ensuring that the UTS is fully extended during the passive mobilization of the operated finger, without making use of the elasticity of the UTS. The UTS may only be cut to the patient's requirements after this extension has been verified.

N.B. We recommend fixing the distal point well because its disinsertion could cause a proximal migration of the UTS. Mounted needles with round points should be used exclusively for this procedure as the use of a needle with a triangular point could cause cuts in the silicone which could then spread.

TO OPEN A PACKAGED STERILE PRODUCT

1. Remove outer packaging under clean, aseptic and sterile conditions and allow the sealed inside packaging to fall gently into place.
2. Attach the registration section of the patient's internal label to the patient's chart.
3. Open inside packaging under aseptic conditions.

EXAMINATION AND HANDLING OF THE PRODUCT

1. The product must be examined visually to detect any signs of particular contamination or other damage.
2. The implant must remain submerged in sterile water or in normal saline solution before its implantation to prevent any contact with external airborne contaminants.

CONTAMINATION

AREX® cannot be held responsible for products contaminated by other substances after the product has left the company.

Surface contaminants (talc, dust, bandages, cotton, wool, oils) on the implants can cause a reaction because they are foreign bodies.

Handle with care wearing surgical gloves (rinsed to remove the talc) and follow directions for aseptic handling.

LIMITED WARRANTY, LIMITED LIABILITY AND RENUNCIATION OF ADDITIONAL GUARANTEES

The manufacturer guarantees that all necessary precautions have been taken during the manufacture of these products and undertakes to replace any product which, following an investigation by the manufacturer, is shown to have been defective at the moment of despatch. The surgeon is entirely responsible for the selection of patients, the surgical procedure, the treatment and post-surgical stress as well as the handling of the product. AREX® has no control over the conditions of use and cannot be liable for the good or bad effect of the product following its use, nor can AREX® be held liable for losses or damages resulting directly or indirectly from the use of this device.

Packaging material disposal should be done according to the law.

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