FR

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EN

PRODUCTS AND USE

- **ALICE pliers** were designed to remove/push and cut K-wires 0.8 mm to 1.2 mm in diameter. **DO NOT USE IT AS A WIRE BENDER.**

- **MANOTTE 1.6TB** pliers can bend and cut K-wires 1.0 mm to 1.6 mm in diameter. They are also used to crimp the LIGAMENTOTAXOR K-wire protection tips. The MANOTTE 1.6TB is designed to install the LIGAMENTOTAXOR.

- **HK2 pliers** can bend and cut 1.2 mm, 1.5 mm and 1.8 mm diameter K-wires. They also crimp the corresponding connectors. HK2 pliers are specifically designed to install the HK2 system.

WARNING

- Do not use the instrument after the sterilization cycle when it is still warm (risk of irreversible jamming).
- Do not allow an instrument to retain humidity on a laser marking. If a brown stain appears, rub hard with a soft cloth.
- Do not store instruments near potentially corrosive substances, particularly chlorine.
- Failure to follow good practices for cleaning and sterilizing surgical instruments may lead to infections or device malfunction.
- The use of a device for a task other than that for which it was designed may lead to premature wear of the device.
- Do not use instruments if they do not perform the intended function correctly or if they are damaged.

LIMITATIONS ON REPROCESSING

The instruments can be reused. However, the operating lifetime of the instruments is limited. The operating lifetime in clinical practice depends on individual use and processing conditions specific to the hospital facility. The end of the operating lifetime is determined by the deterioration and damage related to device use. The devices can be reused when the cleaning and sterilizing conditions recommended by ARES[®] are followed. If one suspects contamination by unconventional transmissible agents or when being used in the context of a high-risk surgical procedure, the device(s) cannot be reused and must be destroyed according to the applicable national regulations.

INSTRUCTIONS

Initial processing at point of use

Always remove any oily residue on items before cleaning.

Containment and transport

Store in a safe, dry place. Transport in a closed container to the reprocessing facility to avoid any damage or contamination to the environment.

Preparation before cleaning

The pliers must always be open before cleaning. It is recommended to clean new instruments before the first sterilization and after each use. The device must be cleaned immediately after use to prevent the drying or fixation of soil.

Immediately after the surgical procedure, the pliers are opened and plunged into a solution containing a cleaning/disinfecting product at a temperature of 45°C for at least 15 minutes.

Do not use cleaning/disinfecting products that fix proteins (i.e. aldehydes). Follow the product manufacturer's instructions on the concentration and duration of action for the cleaning/disinfecting product.

During the validation of the cleaning process, ARES[®] validated the following scenario:

1. The pliers are placed in tap water at room temperature for 15 minutes.
2. The pliers are then brushed with a soft nylon brush (MED 100.33 INSITUMED) under running tap water until all visible residue were removed on the samples tested.
3. The pliers' joints are rinsed with a water spray gun for 15 seconds.

Automated cleaning

Instruments made out of stainless steel, such as ALICE, MANOTTE 1.6TB and HK2 pliers, are designed to withstand intensive cleaning in sodium hydroxide solutions, notably used for elimination of TSE agents. Open the pliers before placing them in the washer/disinfectant to ensure all parts of the pliers will be cleaned.

Steps:

1. **Precleaning**
Tap water: 16°C ± 2°C — 3 min
Products used by ARES[®]: No products used.
2. **Drip dry**
3. **Cleaning**
Tap water: 60°C — 5 min
Products used by ARES[®]: Deconex prozyme alka (Borer). Concentration: 0.5%
4. **Drip dry**
5. **Intermediate rinse**
De-ionized water: 45°C — 1 min
Products used by ARES[®]: No product used.
6. **Drip dry**
7. **Final rinse**
De-ionized water: 20°C ± 2°C - 0.5 ml/l — 1 min
Products used by ARES[®]: No product used.
8. **Drip dry**
9. **Thermal disinfection**
Demineralized water: 90°C — 5 min
Products used by ARES[®]: No product used.
10. **Drying**
No water used: 90°C — 15 min
Products used by ARES[®]: No product used.

Two automated cleaning cycles can be done before sterilization when needed. After cleaning, always make sure the device(s) is clean, dry, functions properly and not damaged.

Manual cleaning

- Steps:
- 1. **Precleaning**
Tap water: 16°C ±2°C — 5 min
Products used by AREX®: No product used.
 - 2. **Precleaning**
Ultrasonic bath: 20°C — 5 min
Produits utilisés par AREX®: Deconex Prozyme Active (Borer). Concentration: 0.2%
 - 3. **Rinsing**
Tap water: 16°C ±2°C — 20sec
Water spray gun with static pressure of about 2 bars.

Two manual cleaning cycles can be done before sterilization when needed.
After cleaning, always make sure the device(s) is clean, dry, functions properly and not damaged.

- Drying**
- **Automated drying**
Dry the pliers using the washer/disinfector's drying cycle as described above. If needed, additional manual drying can be done with a lint-free towel.
 - **Manual drying**
Dry the pliers with a lint-free towel.

MAINTENANCE
After cleaning, apply a small amount of surgical instrument lubricant on the joints for proper functioning.

- INSPECTIONS AND FUNCTIONAL TESTS**
- After every cleaning, and drying, check the following on the product: dry state, cleanliness, proper functioning, and absence of damage.
 - Check that joints move easily, without excessive play. Check the function of locking mechanisms (hook, clamping, etc.).
 - Inspect visually for signs of damage and wear. Cutting edges must be continuous and unbroken.
 - Check the pliers' gripper axis to detect any sign of distortion.

WRAPPING
Standard individual packaging compatible with steam sterilization must be used (i.e. sealed pouches, box with filter and valve or crepe paper). Make sure that the pouch is sufficiently large to house the instrument without stressing the bag's weld.
During the validation of the sterilization process, AREX® used double sterilization envelopes (STERIKING® flat roll manufactured by WIPAK). The pouches were sealed by a HAWO brand heat sealer.

STERILIZATION
Open the pliers before placing them in the sterilizer. The ALICE, MANOTTE 1.6TB and HK2 pliers must be sterilized in a vacuum autoclave using one of the following protocols:

- Protocol A**
- 3 pre-vacuum phases
 - Sterilization temperature: 132°C
 - Holding time: 4 min
 - Drying time: 18 min

- Protocol B**
- 3 pre-vacuum phases
 - Sterilization temp.: 134°C
 - Holding time: 3 min
 - Drying time: 18 min

The autoclave manufacturer's instructions must be followed.
Make sure that all of the device's inner and outer surfaces are freely accessible during the sterilization.
Only cleaning and sterilization processes defined in this leaflet have been validated; however, the healthcare unit may use their own standard disinfection and sterilization protocol, provided the said protocol has been validated for the disinfection and the sterilization of surgical instruments.

STORAGE CONDITIONS
ALICE, MANOTTE 1.6TB and HK2 pliers must be stored in a clean and dry environment, at room temperature. Do not store the instruments near products that may have a corrosive action (chlorine).

ADDITIONAL INSTRUCTIONS
ALICE, MANOTTE 1.6TB and HK2 pliers are **supplied NON-STERILE**.
They can only be used by qualified personnel and must only be used for the type of surgery for which they were designed. Any breach of this rule voids our warranty: 5 years for instruments, 1 year for parts in use such as jaws.
If for any reason an instrument needs to be **returned** to us, it must be first **decontaminated**. AREX® reserves the right to request a proof of decontamination.

- PRECAUTIONS FOR USE**
- All instruments covered by this instruction for use must be used according to this IFU. All information related to the use of these instruments is available on AREX® web site: www.arex.fr
 - Read all sections of this leaflet before use.
 - If the instruments are broken or damaged, they must be discarded according to applicable national regulations.
 - It is recommended to clean the instruments immediately after use.
 - Before every use, inspect the device(s) to ensure its condition and proper function.

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