CHIROBLOC

main malléable / malleable hand
Formbare Hand / apoyamano maleable
mano malleabile / mão maleável
vormbare hand

CHIROBLOC Ballast

lest pour main malléable / ballast for malleable hand
Ballast für formbare Hand / lastre para el apoyamano
maleable / zavorra per mano malleabile
estabilizador para mão maleável
ballast voor vormbare hand

RING

élastique de fixation de doigt pour CHIROBLOC
Elastic finger fixation for CHIROBLOC
Gummi zur Fingerfixierung für CHIROBLOC
elástico de fijación de dedos para CHIROBLOC
elastico di fissaggio del dito per CHIROBLOC
elástico de fixação de dedo para a CHIROBLOC
elastische vingerfixatie voor CHIROBLOC

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DEVICE AND INDICATION FOR USE

• **CHIROBLOC [CBP, CBM and CHIRO3]**

CHIROBLOC, a reusable malleable silicone hand, is designed to hold the patient's hand in the desired position during an orthopedic surgery. CBM is intended to be used on adult size hands, CBP is intended to be used on children size hands, CHIRO3 can be used on all size of hands.

CHIROBLOC is compatible with MRI (*Magnetic Resonance Imaging*).

• **CHIROBLOC Ballast [CBL]**

CHIROBLOC Ballast is a reusable stainless-steel accessory which allows to ballast the assembly in order to stabilise CHIROBLOC in the desired position.

• **RING**

RING is a reusable silicon component of CHIRO3 designed to be used as an elastic finger fixation. RING can be also used as an accessory with all CHIROBLOC range for the same function.

These devices have been designed to facilitate hand surgery. These devices which are supplied non-sterile must undergo an automatic cleaning step followed by a steam sterilization step before use.

WARNINGS

- Do not use the instrument when it comes out of the sterilization cycle and is still warm : there is a risk of premature wear.
- CHIROBLOC must be completely unfolded and well flattened before cleaning and sterilization otherwise its lifetime will be divided by 2.
- Do not let humidity remain on the CHIROBLOC Ballast's laser marking. Otherwise, if a brownish spot appears, rub it away vigorously with a soft cloth.
- Do not store near corrosive products (*ex: chlorine*).
- Do not use the instruments if they are damaged. Example: If a defect is found in the CHIROBLOC's plugs, it is non-compliant and no longer usable.
- Do not clamp RING with a metal instrument, this can make small cuts in the silicone that may cause failure later on.
- Failure to follow good practices for the cleaning and sterilization of surgical instruments may lead to infections or device malfunction.
- Do not use CHIROBLOC and its accessories for any use other than those defined by the manufacturer.

LIMITATIONS ON REPROCESSING

The instruments can be reused, however, their lifespan is limited. The lifespan of CHIROBLOC and its accessories depends on the individual use (*wear, damage*).

The devices can be reused if the conditions of cleaning and sterilization provided by AREX® are followed. The mechanical strength of these devices has been confirmed by AREX® after the equivalent of 225 cleaning/sterilization cycles.

In order to avoid damaging of CHIROBLOC and its accessories, it is forbidden to use solutions containing: chlorine, soda and aromatic derivatives compounds (*Ether, Chloroform, Benzene, Xylene etc.*) or any other abrasive product when the cleaning and the decontamination these devices. In the case of suspected contamination by non-conventional

transmissible agents (*NTCA: prion disease*), the devices cannot be reused and must be destroyed according to current regulations.

INSTRUCTIONS

1. Putting the new product into circulation (*out of the packaging*):

Remove the devices from their packaging.

It is recommended to clean the new device before the first sterilisation according to point 5.

2. Initial processing at point of use (*after surgical use*):

Cleaning must be carried out immediately after each use of the medical device to prevent corrosion and dirt adhesion. Use a cloth to remove any excess dirt.

3. Containment and transport

The mode of transport must be secured, no additional requirements.

4. Preparation before automatic cleaning

It is recommended to clean new instruments before their initial sterilization and immediately after each use in order to avoid dirt adhesion.

- Unfold CHIROBLOC [*CBP/CBM/CHIRO3*]
- Separate the CHIROBLOC Ballast from CHIROBLOC
- Separate the strap from the CHIRO3
- Separate RING

CHIROBLOC [*CBP/CBM/CHIRO3*], CHIROBLOC Ballast and RING will be cleaned like surgical instruments.

5. Automated cleaning

The used detergent is an enzyme-based or alkaline based.

Use a cleaning program similar to this one described below:

1. Pre-cleaning
Tap water / Cold water / 2 min.
2. Drying
3. Cleaning
Tap water + Medizyme enzymatic ou Mediclean forte, alkaline (Neodisher) 0.5% / 55°C / 5 min.*
4. Drying
5. Intermediate rinsing
Demineralized water / Cold water / 3 min.
6. Drying
7. Final rinse
Demineralized water / Cold water / 2 min.
8. Drying

** Follow the detergent manufacturer's instruction.*

Note: A high temperature disinfection step can take place after the final rinse to prevent the fixation of any remaining proteins.

It is recommended to use a washer/disinfector conforms to ISO 15883.

2 automatic cleaning cycles can be performed prior to sterilization if necessary.

After cleaning, always check that the devices are dry, clean, in good working order and absence of deterioration.

6. Automated drying

Dry CHIROBLOC, RING et CHIROBLOC Ballast according to the drying cycle of the washer. If necessary, additional manual drying can be performed using lint free towel.

7. Maintenance

Non-applicable.

8. Controls and functional tests

After each cleaning and drying, check the dryness, cleanliness, proper functioning and absence of deterioration of the devices.

Before each use, carry out the following tests in order to verify the conformity of the device:

- **CHIROBLOC**

Visual inspection of the device to identify any damage, wear and/or signs of degradation: hole, crack, air pocket in the hand, defect on CHIROBLOC's plugs, broken finger under the silicon cover.

- **CBL**

Visual inspection of the device to identify any damage wear and/or signs of degradation : corrosion

- **RING**

Visual inspection of the device to identify any damage wear and/or signs of degradation: crack

9. Packaging

Standard individual packaging compatible with steam sterilization method must be used (*i.e. sealed bags*). Make sure the pouch is enough large to hold the unfolded instrument without placing stress on the pouch's seal. The CHIROBLOC Ballast is always packaged separately.

10. Sterilization

CHIROBLOC must be unfold and flattened before sterilization. CHIROBLOC Ballast, CHIROBLOC et RING are sterilized in a vacuum autoclave according to one of the following protocols:

Protocol A

- 3 pre-vacuum phases
- Sterilization temperature: 132°C
- Holding time: 4 minutes
- Drying time: 20 minutes

Protocol B

- 3 pre-vacuum phases
- Sterilization temperature: 134°C
- Holding time: 3 minutes
- Drying time: 20 minutes

The instructions of the autoclave's manufacturer must be followed.

Ensure that the surfaces of the device are not covered during sterilisation.

Only the cleaning and sterilization procedures defined in this IFU have been validated. However, the healthcare facility can use its standard disinfection and sterilization protocol, as long as this protocol has been validated for the disinfection and sterilization of surgical instruments (*e.g. sterilization in instrument holders, manual cleaning, other cleaning and sterilization parameters*).

STORAGE CONDITIONS

CHIROBLOC [CBP/CBM/CHIRO3], CHIROBLOC Ballast and RING must be stored in a dry and clean and dry environment, at room temperature. Do not store the instruments near products that may have a corrosive action (*ex: chlorine for CHIROBLOC Ballast*).

DISPOSAL

Any instrument that is broken, damaged or at the end of its life must not be used. It must be destroyed like any medical waste according to current regulations.

ADDITIONAL INSTRUCTIONS

1. General Information

CHIROBLOC [CBM/CBP/CHIRO3], CHIROBLOC Ballast and RING are provided non-sterile.

Devices must be used by surgeon and qualified personnel. Instruments must be used only for the surgical purpose and type for which they have been designed (see indications). Any failure to respect these rules cancels the guarantee:

- CHIROBLOC [CBM/CBP/CHIRO3]: 1 year
- CHIROBLOC Ballast: 5 years
- RING: 1 year

If for any reason, some instrument needs to be returned to us, it must be first decontaminated. AREX® reserves the right to request a decontamination proof.

2. Precautions for use

- CHIROBLOC [CBM/CBP/CHIRO3], CHIROBLOC Ballast and RING must be used in accordance with this Instruction for Use. All the information related to use of these instruments is available on the AREX® website: www.arex.fr
- Read all the sections of this IFU before use.
- The instruments should be cleaned immediately after use.
- Before use, inspect devices to verify its condition and proper functioning.

3. Incident declaration

Any serious incident that has occurred in relation to CHIROBLOC and its accessories should be reported to AREX®, and the competent authority of the Member State in which the user and/or patient is established.

Reminder: A serious incident or risk of serious incident is any incident or risk of incident involving a medical device which led or might have led to death or serious deterioration in state of health of patient, end-user and third party.

These instructions have been validated by the medical device's manufacturer as being appropriate for preparing a medical device for reuse. The department is still responsible for ensuring the sterilization process achieves the intended result by involving the sterilization facility's staff, materials and equipment. This generally requires routine verification and control of the procedures. Similarly, any deviation from the instructions provided should be evaluated to determine the effectiveness and possible adverse consequences.

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