Symboles . Symbols . Symbole
Símbolos . Simboli . Símbolos
Symbolen

Attention
Caution
Achtung
Atención
Attenzione
Atenção
Let op

Consulter les instructions d'utilisation
Refer to the instructions for use
Gebrauchsanweisung beachten
Consulte las instrucciones de uso
Consultare le istruzioni per l'uso
Consultar as instruções de utilização
Raadpleeg de gebruiksaanwijzingen

Numéro de lot
Batch number
Chargenbezeichnung
Número de lote
Numero di lotto
Número de lote
Batchnummer

Numéro de référence
Catalog number
Bestellnummer
Código de referencia
Codice prodotto
Número de referência
Referentienummer

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EN INSTRUCTION LEAFLET For the personal attention of the surgeon

Mini external fixator - HK2®

GENERAL INFORMATION

HK2® is a mini external fixator used for orthopaedic surgical treatment. The box made up of 20 K-wires and 5 double connectors. One HK2® unit is made up of two stainless steel ISO standard 5832-1 K-wires and one aluminium connector. The device is for single use. HK2® system is delivered non sterile in a heat sealed plastic pouch, and have to be cleaned and sterilised. HK2® is a system allowing different fracture treatment, with two diverging K-wires, together linked with a connector. There are three different connectors: 1.8 mm, 1.5 mm and 1.2 mm. Each connector corresponds to two K-wires with the same diameter.

MEDICAL DEVICES ASSOCIATED TO HK2® SYSTEM

HK2® PLIERS should be used to put the HK2® system in place. They serve to bend and cut the K-wires and crimp the connector over two K-wires appropriately positioned. Furthermore, the K-wires may be set into the bone with a drilling motor. The system is left in place for 28 to 60 days (according to indication and location).

INDICATIONS

Patient population: adults and children*.

After required radiological and clinical examinations of the bone, HK2® can be used in the following cases:

- Distal radial fractures with dorsal tilt (generally with 1.8 mm K-wires and corresponding connectors).
- Particularly open fractures or in the event of skin abrasion.
- New fracture on the existing malunion site of a previous fracture.
- Fractures treated with plates in the presence of major osteoporosis and imperfect stability.
- Metacarpals (generally with 1.5 mm K-wires and corresponding connectors) fractures.
- Finger's phalanx (generally with 1.2 mm K-wires and corresponding connectors) fractures.
- Or: any other fracture that can be treated with two diverging K-wires linked together with a connector in the adequate sizes in the upper limb.

*For children, smaller diameter will be used (1.2 mm, 1.5 mm). The choice of the appropriate diameter will be determined according to individual clinical case and morphology.

CONTRAINDICATIONS

The contraindications include but are not limited to:

- Acute or chronic, local or systemic infections; Fever.
- Known allergy or intolerance to the materials.
- Infectious disease; Systemic and metabolic disorders.
- Drug or alcohol dependency.
- Severe osteoporosis, Local bones tumors, severe stiffness of the joints.
- Poor blood supply of the soft tissues, Insufficient bone stock.
- Patient in inappropriate condition, non cooperative or with mental illness.
- Intense physical activity.

The surgeon must warn the patient of the influence of these factors on the success of the intervention, and give advice for minimizing these risks.

OPERATING MODE

The success of the operation is closely related to compliance with indications to the operating technique. We advise the surgeon to attend a mentor training course (with a surgeon already experimented with the use of this technique), before starting the use of HK2®. The present instructions give general guidelines to use HK2® system. The pinning technique used actually depends on the clinical/surgical case treated, and the set of the external fixation system is the sole competence of the orthopaedic surgeon.

INTRA-FOCAL PINNING

See illustration A

- **Bending:**
1 bent intrafocal K-wire
1 straight bracket K-wire
- **Connection:**
Pins are connected at their external parallel section with 1 connector.

We are describing here below the intra focal K-wiring:

1. Insert an intrafocal K-wire as in Kapandji's intrafocal pinning technique.
2. The intrafocal K-wire is bent at the point where it emerges from the skin so that the external tip of the K-wire is perpendicular to the axis of the bone. Use HK2® PLIERS to bend the K-wire. Always open the pliers as wide as possible to facilitate the procedure. The pliers can bend the K-wire up to 90°. K-wires are bent laterally, perpendicular to main axis.
3. After bending, rotate the K-wire to reach the desired position. Use radioscopic control to verify result and take note of alignment with joint surface.
4. Cut the K-wire off approximately 20 mm away from where it emerges from the skin (after bent portion).
5. Position the connector. Do not put the connector too close to the skin. By rotating the connector, the optimum

site can be chosen for insertion of the distal K-wire.

6. Through the same incision as the one used for the first K-wire, insert a juxta-glenoidal distal K-wire parallel to the glenoidal cortex. This K-wire provides stability to the intra focal proximal K-wire, and support for the reduction of fracture.
7. Control by palpation and two X-rays, one at the start and one at the end of the insertion procedure.
8. The two K-wires are cut a second time to adjust them to the same length, as flush as possible to the connector. Take into account connector and pliers' jaws thicknesses.
9. The connector covers the tips of the K-wires. The intra-focal K-wires and distal K-wire are firmly connected after crimping the connector, using HK2® PLIERS.
10. Repeat these operations on each site depending on fracture topography.
11. Control clinically amplitudes and final aspect. Check clinically and with radiography the stability of the assembling.
12. Place an appropriate dressing at K-wire exit, that avoids any contact between the connector and the skin.

EXTRA FOCAL PINNING

See illustration B1

- **Bending:**
2 bent K-wires to form a "bridge" outside the skin.
- **Connection:**
Pins are connected at their external parallel section with 1 connector.

See illustration B2

- **Bending:**
2 bent K-wires to form an "arch" outside the skin.
- **Connection:**
Pins are connected at their external parallel section with 1 or 2 connectors.

We are describing here below the extra focal K-wiring:

1. If required, reduce fracture using the appropriate manoeuvres.
2. Place a distal K-wire (straight or oblique) through the bone (head, neck, or base). This first K-wire crosses at least 1 cortex, and may anchor other cortex. Control this first insertion procedure by palpation and with radioscopy.
3. Place a proximal K-wire (straight or oblique) through the bone (head, neck, or base). This first K-wire crosses at least 1 cortex, and may anchor other cortex. Control this second insertion procedure by palpation and with radioscopy.
4. Check the good position of the 2 K-wires by fluoroscopy and palpation.
5. Use HK2® PLIERS to bend the K-wires. Always open HK2® PLIERS as wide as possible to facilitate the procedure. They can bend the K-wires by up to 90°. K-wires are bent laterally, perpendicular to the pliers. Bent K-wires, at the point where they emerge from the skin, so that the external bent extremities of both K-wires are parallel to each other. Both K-wires must be externalized 1 to 2 cm beyond the skin. Cut the surplus using HK2® pliers cutting function.
6. If necessary, orient the K-wires by rotating them to the desired position.
7. Introduce K-wire tips into connector's holes to join the 2 K-wires together. Before crimping, distraction, neutral, or compression forces may be applied between the two K-wires. Crimp the connector using HK2® PLIERS.
8. After locking the system, cut the parts of the 2 K-wires that are beyond the connector's level to prevent system crowding. Control clinically amplitudes and final aspect. Check clinically and with radiography the stability of the assembling.
9. Place an appropriate dressing at K-wire exit, that avoids any contact between the connector and the skin.

Regularly check the patient:

It is mandatory to make clinical checks, and with x-rays, after the surgical procedure, to evaluate patient evolution, to avoid complication that might lead to a protrusion of the system, which may result in removing the implant.

Post-operative follow-up:

Inpatient pain management on first night. Step 3 analgesics for first 24 hours post-op. Step 2 analgesics for 10 days. Renew dressings until K-wires are extracted. Encourage mobilisation. Non-systematic physiotherapy. No splint. K-wires ambulatory extraction in office, or in theatre, without anaesthesia on day 28 or up to day 60 (according to indication and location). Cut the connector and remove K-wires by traction, or unscrewing the threaded K-wires.

POSSIBLE COMPLICATIONS AND WARNINGS

The compatibility of the HK2® external fixator in the MRI environment has not been evaluated. MRI will not be performed on a wearer of an HK2® external fixator because this type of examination could lead to burns. Use of appropriate surgical procedures and techniques is the sole competence of the medical profession. Some degree of risk is inherent to all surgical procedures and complications are sometimes unavoidable. It is the surgeon's responsibility to explain fully the benefit-to-risk ratio of the procedure under consideration. Potential complications and caveats must be discussed with patients before a decision is taken to proceed with surgery. Patients must be informed of warning signs which should lead them to call their surgeon for remedial action to be taken: system failure (loosening), traction failure, unexpected pain, infection, purulent discharge.

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Complications which may arise as a result of all surgical procedures must be discussed with patients: infection, haematoma, accumulation of fluids, loss of sensitivity, possible allergic reaction to medications used during surgery, or to any component (*chrome and aluminium in particular*), nerve damage, secondary displacement, tendon rupture, patient intolerance to a foreign body, imperfect healing. Question patients thoroughly regarding previous allergies before operating.

SINGLE USE

HK2® are designed for SINGLE USE. The product should be examined visually to detect any particular contamination or other damage. Once used and explanted, they must not be re-used (*K-wires twisted and cut, connector crimped to K-wires*).

INSTRUCTIONS FOR PROCESSING

Preparation before cleaning

HK2® is delivered NON-STERILE. Must be cleaned with appropriate products before use. Remove HK2® from its wrapping and check for any damage. Do not use any damaged product.

Cleaning, disinfection and drying

The device will undergo automated cleaning. Use cleaning program similar or equivalent to the one described below.

Cleaning process:

- 1. Pre-cleaning
Tap water < 45°C / 2 min.
- 2. Cleaning
Neodisher MediClean Forte 0.4%: Heated 55°C / 5 min.*
- 3. Rinse 1
Tap water < 45°C / 2 min.
- 4. Rinse 2
Tap water < 45°C / 2 min.
- 5. Thermal Rinse
Reverse Osmosis Water / 93°C / 5 min.
- 6. Dry
> 50°C / 22 min.

* Recommendations of detergent manufacturer will be followed.
The used detergent is alkaline based solution.

Inspection

After cleaning and drying, check on the product the following elements: dry state, cleanliness, proper functionality, and absence of deterioration.

Packaging

Individual: standard packaging material compatible with the moist-heat sterilization process must be used (*i.e. welded bags, box with filter or valve, ...*). Make sure that the pouch is sufficiently large to house the implant without stressing bag welding.
When validating the sterilization process, ARES® used FDA compliant SPS Medical double pouches.

STERILIZATION

After cleaning, HK2® system must be sterilized before use in pre-vacuum autoclave according to the following protocol.

Sterilization Process

- Cycle Type: Pre-vacuum
- Set Point Temperature: 134°C
- Full Cycle Time: 3 min.
- Drying time: 20 min.
- Cool time (*inside and outside the autoclave chambers*):
Maximum 60 minutes

Make sure that the product of sterilization has access to all inward and outward surfaces. In case of sterilization of several products in a vapour autoclave, make sure that the indicated maximal load authorized is not exceeded.
Only cleaning and sterilization processes defined in this leaflet were validated, however, the relevant Healthcare Unit may use their standard disinfection and sterilization protocol, provided the said protocol is validated for the disinfection and the sterilization of implants.

STORAGE

HK2® must be carefully stored, in a clean place, under normal temperature and moisture conditions. The device must be protected from any corrosive environment.

TRANSPORTATION TO POINT OF USE

According to the validated protocol of the department.

CONTAMINATION

ARES® cannot be held responsible for products contaminated by substances once the product has left its premises. Surface contamination (*talcum powders, dust, dressings, wadding, oil*) on the devices may trigger reactions since they are foreign bodies. Handle with care using surgical gloves (*rinsed out to eliminate talcum powder*) taking care to preserve asepsis at all times.

ELIMINATION

The explanted device must be eliminated according to the national applicable laws.

LIMITED WARRANTY, LIMITATION OF LIABILITY AND WAIVER OF OTHER GUARANTEES

The manufacturer guarantees that all necessary precautions were taken during manufacture of the products and undertakes to replace any product which, following the manufacturer's investigation, proves to have been defective when it was shipped. Patient selection, surgical procedures, treatment, post-operative stress and handling of the system are entirely under surgeon's responsibility. ARES® has no control over conditions of use nor on the good or bad effects of the product's use. ARES® will not be held

responsible for loss or damages as a direct or indirect result of using the system.

INCIDENT DECLARATION

Any serious incident appearance related to HK2® and its accessories must be notified to ARES®, and to the local competent authorities where the user/patient is established.

Reminder: A serious incident denotes any direct or relevant indirect link between the use of a medical device and the death of a patient or user, a «serious public health threat», or the temporary/permanent loss or deterioration of a patient or user's state of health or wellbeing.

Those instructions were validated by the manufacturer of the device allowing to prepare a device in preparation for its reuse. It falls to the competent service to ensure that the cleaning and sterilization process attained the expected result. Those process were installed with the equipments, the materials and the personnel of the sterilisation facility. It generally requires the verification and the routine control of the process in question.

Réf.: IFU_HK2 Ed.01 - 09/2020

Revision date of the instruction for use September 2020.

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