

LIGAMENTOTAXOR REF.: LTX - Kit sold sterile

Dynamic External Fixator

Possible palmar sub-luxation with a third K-Wire,
at the same time, ligamentotaxis change is possible.




PIP joint crush fractures, 3rd & 4th fingers

Male - 54 years —

By Philippe Pélissier, MD

Bordeaux University Hospital - France



www.arex.fr/QRc.php?c=16



Quick video to present LIGAMENTOTAXOR

www.arex.fr/QRc.php?c=1

MANOTTE 1.6TB REF.: MAN16TB - sold non sterile

The instrument for LIGAMENTOTAXOR technique.
To cut and bend 1 mm to 1.6 mm diameter Kirschner wires.



Detailed video to present MANOTTE 1.6TB pliers

Bend / Cut / Crimp

www.arex.fr/QRc.php?c=8




AREX®: an innovative company dedicated to hand surgery. A vast range of specialised instruments, specific implants and surgical techniques.

AREX® INSTRUMENTS - reusable, sold non sterile


ALICE pliers REF.: ALICE

- A wire puller to simplify a difficult technique for KW 0.8mm to 2 mm.
- A wire cutter: 1.2 mm ø max



ForgetMeNot REF.: FMN Y-2B (yellow) or FMN B-2B (blue)

- Simple answer to a real issue: digit tourniquet ForgetMeNot
- Available in yellow or blue
 - Sold in packs of 10




CBE - CARPAL BONE EXTRACTOR REF.: CBE

- It is not a cork-screw but a joystick helping moving a carpal bone you want to remove.



MATBOX range REF.: MBX / MBXII

- A sterilization container, to protect and store fragile surgical instruments.
- Usable dimensions:
MATBOX: 200x200x20mm / MATBOXII: 100x200x20mm



DISTRIBUTION PRODUCTS

MICROSURGERY INSTRUMENTS

- Several models available.
- More information on our website: www.alex.fr
- Sold non sterile.



MICRO VASCULAR CLAMP

- Single use.
- Sold sterile.

Distributed in France by AREX®, for other countries please contact BIOVER®.

TKFV-1 15 g/mm ²	TKSV-1 20 g/mm ²	TKMV-1 30 g/mm ²	TKLV-1 40 g/mm ²
TKFV-2 15 g/mm ²	TKSV-2 20 g/mm ²	TKMV-2 30 g/mm ²	TKLV-2 40 g/mm ²

TKS-1 40 g/mm ²	TKM-1 60 g/mm ²	TKL-1 120 g/mm ²
TKS-2 40 g/mm ²	TKM-2 60 g/mm ²	TKL-2 120 g/mm ²

SCRU[™]2[©] and DIP Fusion

Low Profile and Standard with identical length have exactly the same pitch, thus Standard SCRU[™]2[©] may act as a « spare screw » for LP SCRU[™]2[©].

SCRU[™]2[©] **Low Profile** REF.: P2PxxLP - sold non sterile
2.3 mm \varnothing top of the threads below bed nail
guided over a 1 mm \varnothing KW



SCRU[™]2[©] **Standard** REF.: P2Pxx - sold non sterile
2.5 mm \varnothing top of the threads below bed nail
guided over a 1 mm \varnothing KW




Quick video to present surgical technique



www.alex.fr/QRc.php?c=15



Michel Merle, MD



FIRST THREAD
• 2.5 mm \varnothing top of the thread



SECOND THREAD
• 2.3 mm \varnothing top of the thread

SDNS REF.: SDNS - sold non sterile

The instrument for AREX SCRU[™]2[©] handling.



Universal Tendon Spacer REF.: TCTU1 - sold sterile


UTS & metacarpal hand replantation

35 year old male, transmetacarpal amputation of the hand, circular saw accident.

Flexor tendon avulsion injuries in amputated fragment.

by Philippe Pélissier, MD - Bordeaux University Hospital, France

Replantation and 1st STEP of two-stage multiple flexor tendon grafting




Metacarpal shortening to enable replantation, combined with first stage flexor tendon reconstruction using Hunter's technique, UNIVERSAL TENDON SPACERS and preparation of flexor tendon grafts by the Paneva-Holevitch technique (suture of flexor digitorum profundus [FDP] over flexor digitorum superficialis [FDS]).

Replantation and 2nd STEP of two-stage multiple flexor tendon grafting

The second stage combines flexor tendon grafting with extensor system tenolysis.

8 months after 1st STEP



6 months after 2nd STEP



12 year follow up video



12 year follow up video

www.arex.fr/QRc.php?c=17



HK2 - Mini External Fixator according to Metahus technique




ISELIN TECHNIQUE MODIFIED BY PH. LIVERNEAUX

Complete technique video

www.arex.fr/QRc.php?c=20

Case from Philippe Liverneaux, MD PhD - Strasbourg University Hospital - France



Video 6 months post-op



www.arex.fr/QRc.php?c=18

Case from Christian Thomas, MD - Centre de la Main - Toulon, France



Video 40 days post-op



www.arex.fr/QRc.php?c=29


- One HK2 unit** - sold non sterile
- 2 K-wires, same diameter (1.2 mm or 1.5 mm or 1.8 mm)
- 1 double connector to join the 2 K-wires together

One box - 10 HK2 units
 Three different references
 REF.: HK212 or HK215 or HK218



HK2 pliers REF.: HK2P - sold non sterile

The instrument for HK2 technique.
To cut and bend 1.2/1.5/1.8mm Kirschner wires.



Bend KW without any stress to the bone.

Cut pin end clean & burr-free.

Crimp without modifying KW respective position within connector.

CHIROBLOC3 - Leadless & Radiolucent

Needle Fasciotomy for Dupuytren Contracture

Clinical case from Stéphanie Gouzou MD, Strasbourg University Hospital - France



Quick video to present surgical technique

www.arex.fr/QRc.php?c=21



CHIROBLOC 3 is sold with 5 RING

RING is used for hand positioning with all CHIROBLOC models.

3 models of CHIROBLOC with different ergonomies

Reusable, sold non sterile



CHIROBLOC 3

REF.: CHIRO3



CHIROBLOC Large

REF.: CBM



CHIROBLOC Small

REF.: CBP

ACCESSORIES

Reusable, sold non sterile



BALLAST

REF.: CBL



RING

REF.: RING

Sold in 10 units pack

AREX[®]

3 allée du clos tonnerre 91120 PALAISEAU - FRANCE

Tel. +33 [0] 169 412 212

Fax +33 [0] 169 411 110

info@arex.fr

www.arex.fr

Asia Pacific: **AREX[®]ASIA**

Tel. +603 5121 9823 - arexasia@gmail.com