

Anywhere you can use a K wire you can use a **chiroKlip**™

- Bend K wire as usual, apply self-adhesive chiroKlip™ over sharp end.
- Prevents injury.
- Stops snagging of K wire.
- Makes removal of dressing easier - chiroKlip™ will not catch in the dressing.
- For fingers, toes - or anywhere external.
- Allows for easy removal of K wire, no pliers needed.
- Patent applied for: GB 0123670.2
- Not for implantation. For external use only.


CE Manufacturer:
SYNERCEPTION Limited
Unit 1, 80 High Street
Winchester
Hampshire
SO23 9AT - United King

www.arex.fr

AREX® - 91120 PALAISEAU - FRANCE — T. +33 [0] 169 412 212 _ F. +33 [0] 169 411 110 _ info@arex.fr

Design: e-monceaux.fr / 096 64.04 / © AREX® - 2019

Anywhere you can use a K wire you can use a **chiroKlip**™



- Protects Kwire
- Easy to attach
- Aids easy removal
- Prevents needle stick injuries

• Twenty five handy packs of two

• No special instruments required

• Sterilized for immediate use