


LIGAMENTOTAXOR

Clinical case, Prof. Ph. Pélissier / Bordeaux CHU, France

Male - 54 years

Crush fracture of proximal interphalangeal (PIP) joints, 3rd and 4th fingers.



On 3rd finger, articular traction with LIGAMENTOTAXOR to reduce fracture.

On 4th finger, initial posterior displacement at the base of P2 is reduced by a supplementary transverse pin supported by lateral bars.

Post operative

Rehabilitation protocol:

D00 - D10 Articular traction, no mobilisation.

D10 - D25 Passive mobilisation.

D25 - D40 Active mobilisation.

D40 Apparatus is removed.

Results after 4 months:



Video at 18 months



<http://www.arex.fr/QRc.php?c=16>